



EMTrack™ Release Notes

These notes describe highlights of the EMTrack 4.3 release that is scheduled for Wednesday, April 24, 2019. Additional details about these and other features are available through your solution. For more information, please contact your Client Success Manager (CSM) or Juvare Support at 877-771-0911 or support@juvare.com.

Enhanced Incoming Patient Notification Audible Alerts

Enhanced audible alert options allow Regional Administrators to choose the sound associated with incoming patient notifications (IPN). Through *Settings*, administrators can also suppress audible alerts for users that are viewing the *Transports* page. While on this page, users can temporarily mute or unmute IPN alerts for their current session.

Removed “Patient Delivered” Notifications

To reduce the number of notifications received, EMTrack will no longer send “patient delivered” notifications through email, text, and/or pager channels.