



EMTrack® 5.2.2 Release Notes

These notes describe enhancements introduced in the EMTrack 5.2.2 release scheduled for Wednesday, March 3, 2021. More information is available through your solution or by contacting the Juvare Support Center at 877-771-0911 or support@juvare.com.

Tracking Number

Administrators can now set a minimum length for the patient identification or tracking number. This improvement allows you to establish and enforce basic numbering standards.

Ethnic Group

For the Ethnic Group question on patient forms, an option has been added that allows you to identify a person's group as "unknown".

Resolutions

Resolved a minor issue associated with bulk syncing of information with third parties.