

# Contaminated Patient List Report

[Home](#) > [EMTrack](#) > [EMTrack Reports](#) > Contaminated Patient List Report



The Contaminated Patient List Report displays clients who have been reported as contaminated. That is, clients are included if a user selected Yes as the value for the Contaminated field in the client record.

Setting up this report includes selecting a date range and report format.

## To generate this report

1. Click the **Reports** tab.
2. Click **Contaminated Patient List Report**. The *Contaminated Patient List Report* window opens.
3. Select or enter the date and time range (**Start time**, **End time**).
4. Select the **Format**.
5. Click **Run Report**. The *Generating report* window opens.
6. Do you want to generate another report, click **Yes** or **No**.
7. Open or save the report.